

MARS Bar BIBO Tubing Set

Certificate of Compliance



Issue date: 05/23/2025

Product Name: MARS Bar BIBO Tubing Set

Catalog: BRB02

Lot. No: SC25001

Use-by Date: One Day of Use for One Experiment

Storage: All components should be stored at room temperature

Expiration Date: 05/20/2026

Required Testing	Specification	Applied Standards	Result
Cytotoxicity Test	No Cytotoxicity	ISO 10993-5:2009	PASS (No Cytotoxicity to L-929 cells)
Transportation Study Test	Stress Test	ASTM D4169 DC13	PASS
Endotoxin Test	Bacteria Endotoxin <2.5 EU/device		PASS
1-Year Sterility Test	Accelerate Aging Test, Visual Inspection, Vacuum Leak Test, Bust Test,	ISO 11607	PASS
1-Year Accelerated Aging Test	No microbial growth	ISO 11737-2:2019	PASS
PTT Test	No difference in PTT between the sample and the blank control or the negative control		PASS
In Vitro Hemolytic Test	Nonhemolytic		PASS
Acute Systemic Toxicity	No Acute Systemic Toxicity		PASS
In Vitro Test		ISO 10993	PASS

Intracutaneous Reactivity Test	No Adverse Skin Reaction		PASS
Barrier Seal Test	Packaging Integrity	ISO 11607	PASS
Pyrogen Test	Absence of pyrogens		PASS
Sterility Test	Sterile No microbial growth	ISO11737-2:2019	PASS
Animal Content	No materials of animals or insect origin		PASS

This product lot has passed Applied Cell's QC testing. For details on QC testing contact us at info@appliedcells.com

Applied Cells Inc. guarantees the performance of this product until the expiration date. The guarantee is dependent upon proper storage and handling as instructed on our Product Data Sheets. Every lot of products is quality tested against a "gold standard" reference lot. A new lot is only released based on our defined QC specifications to ensure lot to lot reproducibility and reliability. Applied Cells Inc. guarantees the stability and performance of all our products shipped at room temperature.

Prepared By

Applied Cells Inc.
3350 Scott Blvd, Bldg. 6
Santa Clara, CA 95054
Tel: 1-800-960-3004

General Disclaimer

For research use only. Not for diagnostic use. Applied Cells Inc. will not be held responsible for patent infringement or other violations that may occur with the use of our products.